

FFP2 NR PREMIUM RESPIRATOR INSTRUCTIONS FOR USE

Name of the personal protective equipment:

FFP2 NR Premium Respirator

Supplier:

Cool agency, s.r.o, Koněvova 2660/141, 130 00 Prague 3 – Žižkov, Czech Republic

Producer's name and address:

Intended purpose of use: The intended purpose of use of the respirator is to protect its user by filtering the inhaled air and decrease the transmission of COVID-19 infection between persons. The filtration efficiency is $\geq 95\%$. It is a single-use product.

Efficiency: FFP2 NR class.

Special attention should be paid to warning signs.

Application procedure:

Put the respirator on only with clean hands. All parts of the respirator should be checked for damage before use. Use both hands.

1. Unpack the respirator and take the side bands.
2. Place the respirator on your face so that it covers your nose and mouth and ends down under your chin.
3. Fix the bands behind your ears.
4. Check and shape the nose clip with your fingers.
5. Check the respirator for a tight fit around the perimeter.

Illustration picture of the application procedure:



Checking the tightness and correct placement:

1. Cover the front part of the respirator with both hands. Be careful not to change the respirator position.
2. Exhale sharply.
3. If air leaks around your nose, adjust the shape of the nose clip. Check the tight fit again as described above.
4. If air leaks through the edges of the half mask, adjust the side bands around your ears until you eliminate the air flow. Check the tight fit as described above.
5. If the respirator cannot be fitted properly, replace it and do not continue using the original.
6. If air does not leak, you can use the respirator.

Instructions for storage:

Protect from sunshine and moisture.

FFP2NR Premium Respirator meets the requirements of the European Standard EN149:2001 + A1:2009.

Warning:

- Always make sure that the product is complete, properly put on and worn during the whole time of exposition.
- This product protects from gases/vapours.
- Do not use this product as a protection against atmospheric contaminants, the concentration of which is no longer known or which are directly harmful to the health or life.
- Do not use this product if you have a beard or any other facial hair which could prevent a tight fit between your face and the protective product.
- This product does not eliminate the risk of contracting a disease or infection.
- Please report any serious incident related to this product to Cool Agency, s.r.o. and a respective local supervisory authority.
- Destroy and replace the FFP2 NR Premium Respirator, if it gets contaminated with blood or other infectious contaminants, if it is

damaged, or if breathing becomes difficult. After the end of its use ensure its disposal in accordance with the requirements of the directive on the disposal of hazardous infectious waste.

- This product is for single use only. Its repeated use can lead to cross-contamination, compromise the user safety and impair the overall performance of the product. Contaminated products must be disposed of as hazardous waste in accordance with national regulations. Do not modify, clean or repair this product.
- This product is not fire-resistant. Do not use it in direct contact with open fire.
- Do not use it in atmosphere containing less than 19.5 % of oxygen (definition of Cool Agency, s.r.o; some countries can have their own definition of an oxygen deficient atmosphere; in case of doubt, seek professional help).

Failure to meet the instructions above or unprofessional handling of the respiratory protective product during the stay of the user in a contaminated environment can have adverse health effects and cause serious illness.

The month and year of production and expiration of the product are given on the package.

Date of the last revision of the Instructions for Use:

2020-09-22.